

**In the Specification:**

On page 19, line 11, please add the underlined reference:

Cell pellets were lysed in 0.4 ml/tube of Schagger-Von Jagow (SVJ) electrophoresis system treatment solution and treated at 95° C for 5 min. Lysates (100 µl each) were separated by SDS-PAGE run according to Schagger-Von Jagow (see Schagger, H., and Von Jagow, G., Analytical Biochemistry 166:368-379, 1987) along with FITC-Insulin standards, and treated aliquots of the cell-unbound incubation mixture.